 Lighthouse Info Systems Pvt. Ltd. <small>SINCE 1987</small>	Document Name	Object Help Document	
	Menu Code	01.501.02	
	Process Name	Item test Master	
Release Version	1.0	Document No	LISPL/OHD/PP/501/002
Release Date	02-Jan-2016	Pages	Page - 3 - of 4

Cancel	To exit from current working form.
D	To delete the Line Item Record.

Work Step

Step 1 Press 'New' Button.

Step 2 Select 'Item Category' from the LOV (Press F9 for LOV).

Step 3 Select Item (Press F9 for LOV).

Enter Data in 'Item test rate details' block.

Step 4 Select 'Test Code' (Press F9 for LOV).

Step 5 Enter Min Limit.

Step 6 Enter Max Limit.

Step 7 Select tolerances Basis from drop down list.

Step 8 Enter the user seq.

Optional Data Fields

Select Process Code (if applicable)

Select Test item code (Press F9 for LOV).

Select Third Party details from drop down list.


Select Single QC Level from drop down list.

Enter sampling criteria that is going to be used.

Enter Test Method

Input Field Level Information

FIELDS	DESCRIPTION	R/O/C
Item Category	To Select Item category from LOV.	R
Process	Select the Process Code if it mapped and applicable, otherwise skip.	O
Item	Select Item by pressing F9 button for LOV, it will populate the list of items from the selected item category.	R
Test Code	Select the Test Code by pressing F9 which would be applicable for the selected item.	R
Test Detail	Displays the Test Details of Test Code	D

 Lighthouse Info Systems Pvt. Ltd. <small>SINCE 1987</small>	Document Name	Object Help Document	
	Menu Code	01.501.02	
	Process Name	Item test Master	
Release Version	1.0	Document No	LISPL/OHD/PP/501/002
Release Date	02-Jan-2016	Pages	Page - 4 - of 4

Base Limit	Base limit is used exclusively for the purchase and sales material. It is used to set the deduction limit on the lower side. Specify value in this column when if it fixed all the time irrespective of any supplier.	O
Min Limit	Define the standard minimum limit of test code on which the QC Test Results will be compared with actual result to pass / fail.	O
Max Limit	Define to the standard maximum limit on which the qc result will be compared with the actual result for deduction and passing the QC test.	O
Tolerance basis	Select the Tolerance basis from the drop down list (Upper / Lower)	O
User seq.	Select the User Sequence and ensure to maintain the User Sequence "1" for the main test code.	R
Test Item Code	User is required to select form the LOV	C
Test Item Name	Gets displayed on the selection of Item code	C
Third Party Testing	Select from the drop down	C
Single QC- Level	Select from the drop down	C

R = Required, O = Optional, C = Conditional D = Display Field